



ReSound ONE

Hear like no other

Experience ReSound ONE



GN Making Life Sound Better

Our most natural
sound quality



Go back

Experience more



Nature has designed the shape of your ear to perfectly capture the sounds around you. And because the shape of your ear is unique to you, the sound you hear is uniquely yours.

Most hearing devices capture sound from behind the ear only. This isn't where nature intended, so sound ends up feeling standardised, unnatural and incomplete.

ReSound ONE™ with M&RIE changes all of this with an entirely new class of full-featured* hearing aid that uses your unique ear shape to collect the whole sound, just as nature intended.

Introducing a truly individualised, more complete sound experience, designed for you.

ReSound ONE. Hear like no other

*Microphone in ear and two standard directional microphones, directionality features, wireless audio streaming



Go back

Experience more



Just like a
fingerprint,
your hearing is
one of a kind



M&RIE
Microphone & Receiver-In-Ear



Go back

Experience more



Introducing the groundbreaking M&RIE

World's first full-featured hearing aid with M&RIE (Microphone & Receiver-In-Ear) hearing aid with two traditional microphones and a third that sits discreetly inside your ear canal.

IN PRIME POSITION

Placing a microphone inside your ear instead of behind it and using your unique ear shape to collect sound give you the truest sense of space and location of sounds.

MADE TO WORK WITH YOU

Immerse yourself in an individualised experience, with your daily sound environments enriched with direction and depth. You and your hearing aid working as one.

INSPIRED BY ORGANIC HEARING

Our sound philosophy gives your brain everything it needs to process sound, making it effortless to tune in and enjoy conversations, in any situation.

ReSound ONE. Hear like no other. Works with you.



Go back

*Microphone in ear and two standard directional microphones, directionality features, wireless audio streaming.

Experience more



Shape your hearing experience

Greater hearing in any environment

With All Access Directionality, ReSound ONE automatically adjusts to your listening situations.

IN CONTROL

You are always in the best position to hear what's important and still monitor the sounds around you comfortably. It lets your brain decide what to listen to, not your hearing aids.

SIMPLY EFFORTLESS

You'll feel at ease and engage effortlessly with others and your changing surroundings, even noisy or challenging listening situations.



[← Go back](#)

[Experience more →](#)



Great 1:1 hearing experience

In the toughest noise environments, activate Ultra Focus. Like an invisible beam of sound detection, it concentrates on the speech of the person in front of you, without isolating you from other sounds, for comfort and safety.

FOCUS IN ON LIFE

Our powerful technology helps you feel confident enough to have conversations in situations that you might have avoided in the past.

ReSound ONE. Hear like no other



Go back

Experience more




Connect to
people and
your world



Go back

Experience more





Unmatched rechargeability and extensive connectivity to adapt to your needs.

RECHARGE AND GO

ReSound ONE has the longest rechargeable battery life for well over a day's worth of worry-free use. Select from two stylish and portable charging cases, to suit your needs and lifestyle.

STREAMING MUSIC, AUDIO AND CALLS

Stream audio directly from your iPhone or Android™ smartphone to ReSound ONE hearing aids and control them with the ReSound Smart 3D app, all via the stable Bluetooth® Low Energy.

A versatile range of wireless accessories give you even more options to connect to the sounds and devices you enjoy every day, like TVs.

ReSound ONE. Hear like no other



 Go back

Experience more 

Designed for ultimate comfort, discretion and performance

Every detail of ReSound ONE was
designed with one purpose – you.




THE FINISHING TOUCHES

The modern, discreet styling sits lower on the ear, with just one touchpoint at the back. ReSound ONE is our most comfortable and discreet hearing aid ever made.

COLOUR STORIES

Select from a timeless range of natural colours to blend in with your skin or hair tone.

ReSound ONE. Hear like no other



Go back

Experience more



At ReSound we believe life is lived best when it's a natural, effortless flow. This means we take anything that holds you back to heart and will always strive to help you connect to a life without limitations. We have devoted the past 150 years to advancing lives through technical innovations that, inspired by our Organic Hearing philosophy, place you at the centre of everything we do, because hearing and hearing care is as individual as a fingerprint.

ReSound is part of the GN Group – pioneering great sound from world-leading ReSound hearing aids to Jabra office headsets and sports headphones. The GN Group was founded in 1869, employs over 5,000 people, and is listed on NASDAQ OMX Copenhagen.

Find out how you can hear like no other with ReSound ONE.

resound.com.

**Worldwide
headquarters**

GN Hearing A/S
Lautrupbjerg 7
DK-2750 Ballerup, Denmark
Tel.: +45 4575 1111
resound.com



Go back