

# Extensive connectivity to adapt to your needs



ReSound ONE adapts to your lifestyle connecting you to the world around you with direct audio streaming from iPhone or Android™ smartphones.



One-tap sound adjustments on your ReSound Smart 3D™ app and a versatile range of wireless accessories give you even more options to connect to the sounds you enjoy every day, like your favourite TV program.

# More connections for better care

ReSound ONE includes access to ReSound Assist services so you can receive hearing care in a safe and convenient way.

## FACE-TO-FACE ONLINE HEARING CARE

ReSound Assist Live makes it possible to have face-to-face fitting sessions from the comfort of your home. Simply answer the incoming call from your hearing care professional on your iPhone or Android™ smartphone and get fine-tuning adjustments made live during your appointment.

## REMOTE FINE-TUNING ANYTIME

ReSound assist remote Fine-tuning also allows you to send requests for additional adjustments anytime through the ReSound Smart 3D app. Updates will be sent for you to download at your convenience - without the need for an office visit.

# Designed for ultimate comfort, discretion and performance



Built for ultimate comfort and adaptability, the modern styling of ReSound ONE has just one touchpoint on the back of the ear and sits lower for greater discretion. It's available in three models.


ReSound ONE comes in a variety of colours to blend with hair and skin tones, so you can wear it with confidence.



## INDUSTRY-LEADING RECHARGEABILITY

ReSound ONE offers the best rechargeable battery life in the industry for well over a day of worry-free use.

For added flexibility, two different charger options serve as stylish, portable and protective cases, ideal for traveling and storage. Choose the Premium Charger with an onboard battery for up to three days of quick charging on the go, without the need for a power outlet, or the smaller Standard Charger that gives you all the power you need in a small elegant design.



\*iPhone® 11 and newer.



ReSound GN

At ReSound we believe life is best when it's a natural, effortless flow, so we take to heart anything that holds you back. This is why we have devoted the past 150 years to making life sound better, connecting you to a life without limitations. Our breakthrough innovations, inspired by our Organic Hearing philosophy, all have one core ingredient – you, because hearing and hearing care should be as individual as a fingerprint.



M102791GB-20.06-Rev.A



# Introducing ReSound ONE

Just like a fingerprint, each person's hearing is one of a kind. This is because nature has designed the shape of your ear to perfectly capture the sounds around you.

ReSound ONE™ with M&RIE (Microphone and Receiver-In-Ear) delivers a new class of hearing aid\* with a third microphone inside your ear, so you can collect the whole sound, just as nature intended.

Add direct audio streaming and industry-leading rechargeable battery life, and ReSound ONE truly keeps you connected to people and technology to become one with your world.

\*Microphone in ear and two standard directional microphones, directionality features, wireless audio streaming

ReSound ONE provides you with a more complete and individualised sound experience with:



### A MORE NATURAL SOUND AND IMPROVED SENSE OF SPACE

ReSound ONE with groundbreaking M&RIE\* is the first hearing aid with three microphones, including one in the ear, to use your ear's natural ability to collect the information your brain needs to make sense of sound.



### EASIER-TO-FOLLOW CONVERSATIONS IN DYNAMIC ENVIRONMENTS

All Access Directionality gives you access to all the sounds around you, putting you in charge of which sounds you want to focus on - not your hearing aids.



### GREAT 1:1 SPEECH UNDERSTANDING

Like an invisible beam of sound detection, Ultra Focus lets you decide when to concentrate on the speech in front of you while minimising background noise - even in the toughest environments.

### Worldwide Headquarters

**GN ReSound A/S**  
Lautrupbjerg 7  
DK-2750 Ballerup  
Denmark  
Tel.: +45 4575 1111  
resound.com

CVR no. 55082715

© 2020 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC.

# ReSound ONE

# Hear like no other

GN Making Life Sound Better